

COMMERCIAL SERIES LIFTS TRI-POINT ANCHORING SYSTEM

Patent Pending

INSTALLATION MANUAL



COMMERCIAL SERIES TRI-POINT ANCHORING SYSTEM INSTALLATION

Patent Pending

PARTS LIST

- ✓ 1 - Anchor Plate with bonding Lug
- ✓ 3 - 5/8 Stainless Steel Bolts
- ✓ 3 - 6" Anchoring Rods



ANCHORING PLATE, ANCHORS & BONDING LUG




5/8's Bolts




ADA GUIDELINES FOR POOL LIFT PLACEMENT


- Pool Lift Location. Pool lifts shall be located where the water level does not exceed 48 inches (1220 mm).
Seat Location. In the raised position, the centerline of the seat shall be located over the deck and 16 inches (405 mm) minimum from the edge of the pool. The deck surface between the centerline of the seat and the pool edge shall have a slope not steeper than 1:48.



- Clear Deck Space. On the side of the seat opposite the water, a clear deck space shall be provided parallel with the seat. The space shall be 36 inches (915 mm) wide minimum and shall extend forward 48 inches (1220 mm) minimum from a line located 12 inches (305 mm) behind the rear edge of the seat. The clear deck space shall have a slope not steeper than 1:48.




- Submerged Depth. The lift shall be designed so that the seat will submerge to a water depth of 18 inches (455 mm) minimum below the stationary water level. To make sure the lift has enough clearance the lift should be installed at a location with 44"-48" of Water depth



RETRO-FIT/COMMERCIAL SERIES



COMMERCIAL SERIES – TRI POINT ANCHORING SYSTEM



Anchor Plate's front edge with must be 17" from the water's edge.

Attach copper wire to the bonding lug and then to your bonding grid

Cold Pin's should be 1/2" rebar 6" to 8" long

NOTE: Make sure your anchor plate is level and square with your pool's edge.



Step 1: Determine the location for installation (Make sure all ADA requirements are met)

Step 2: Make sure that all the parts are present (Refer to the parts list)

- ✓ 1 - Anchor Plate with Bonding Lug
- ✓ 3 - 5/8 Stainless Steel Bolts
- ✓ 3 - 6" Anchors

Step 3: At the determined location for installation, chalk out a 3' X 3' area prior to cutting out the concrete. (Refer to overview, the front end of the Anchor Plate must be 17" from the water's edge and 6" deep). New Construction, the concrete must be 6" deep with a 1/2" reinforcing rod, minimum of an 18" grid.

NOTE: Contact Qualified Contractor for installation!

Step 4: Once the 3" X 3" area is cut out, make sure you drill in the inner wall of all 4 sides of the area (refer to the section detail) for the purpose of cold pinning the new concrete to the old concrete. Use 1/2" rebar, 6" to 8" in length.

Step 5: Attach the 3 – 6" threaded anchor's with the 5/8's bolts (included) to the Anchor Plate and tighten down.



Step 6: Place the Anchor Plate in the area, make sure that the front of the Anchor Plate is 17" from the water's edge. (refer to overview) A key element is to make sure that the Anchor Plate is level and square with the pool so that your Commercial Series Lift sits properly on the deck.

Step 7: Grounding the unit – The bonding lug is already installed on the Anchor plate (furthest point from the pool wall on the Anchor Plate) locate the bonding lug, attach a copper bonding wire from the bonding grid.

NOTE: Contact electrician to make sure you comply with the BONDING CODES.

Step 8: Pour in your concrete, (make sure you cover the holes in the 3 threaded anchor's so you don't get concrete in the threads) minimum of 4,000 psi with reinforcing rod. Also make sure that the tops of the threaded anchors are flush with the deck. Before you install your lift make sure you let the concrete set up for 48 hours.

Step 9: After the concrete has set for 48 hours, now it's time to install your lift. Place your Commercial Series Lift over the anchor holes, screw in the Anchor bolts that are included and tighten down securely. Once you have completed the installation, please refer to User's Manual for a safe operation.